



How to Use Your CalmWand



least calm → most calm



when center light is **ON**,
breathe **IN**



when center light is **OFF**,
breathe **OUT**

EXPERIMENT
try other calming
strategies with your
CalmWand

SAVOR
savor the memory of a
positive experience


REASSURE
put your hand on
your heart, slow
your breathing, and
repeat "i am safe"

BREATHE
breathe in when the
center light is **ON**,
breathe out when
the center light is
OFF

WARM HANDS
repeat "my hands and
legs are heavy and
warm." "my
heartbeat is calm
and regular."

LISTEN
listen to calming
music

VISUALIZE
imagine being in a
beautiful warm
place, or elsewhere
where you feel safe.



GOOD THINGS:
lights change color in a **clockwise** direction
vibration when the lights change
If you notice these things, it shows that you
are heading in the right direction!

How to Use Your CalmWand



Turn it on

We recommend two daily sessions of 10 minutes each.



Have hand sensor facing the palm

Studies show you will sleep better and be more calm with regular use.



Close hand comfortably

With practice, you will learn what it feels like to be physiologically calm and how to get there -- even without the CalmWand!

Index	Temperature Range (°C)	RGB Description
0	70.00 – 72.15	Gray
1	72.15 – 74.31	Dark purple
2	74.31 – 76.46	Red
3	76.46 – 78.62	Orange-red
4	78.62 – 80.77	Orange
5	80.77 – 82.92	Yellow
6	82.92 – 85.08	Yellow-green
7	85.08 – 87.23	Green
8	87.23 – 89.38	Green-cyan
9	89.38 – 91.54	Cyan
10	91.54 – 93.69	Blue
11	93.69 – 95.85	Indigo
12	95.85 – 98.00	Violet
13	> 98.00	White (upper cap)